



TABLE TOP LOW SPEED CENTRIFUGE

SA-C5E

MADE in ITALY



www.sanymed.it



- Microprocessor control, digital display, touch panel, parameters can be edited in running and change to display other running parameter and RCF.
- Brushless frequency motor with simpler construction more reliable performance, longer life and quietly running
 - Automatic lid lock, super speed, over temperature protection and imbalance protection. The centrifuge body is made of high quality steel, safe and reliable.
 - Adapters are available by experiment requirements.
 - 3 tiers protection steel cover, safe and reliable.
 - TUV CE EMC certificate.

CLINICAL CHEMISTRY REAGENTS MANUFACTURER

• TECHNICAL PARAMETERS:

Max Speed	5000 rpm
Max RCF	4420xg
Max Capacity	4x180ml
Time Range	0~99 min
RPM/RCF Covert	Yes
Noise (dB)	≤58
Temperature	Normal
Acc/Dec	10 Kinds
Speed Accuracy	±20r/min
Temperature Accuracy	/
Voltage (V/Hz)	AC 220V/110V 50HZ/60HZ
Size (WxDxHmm)	540x370x280 mm
Net Weight (kg)	28Kg
Certificates	CE ISO & Calibration report are available

• MATCHED ROTORS

Order No	Rotor	Max Speed (rpm)	Max Volume (ml)	RCF
5E-1	Swing Rotor	4000	4x180ml	2580
5E-2	Microplate Rotor	4000	2x3x96 well	2020
5E-4	Swing Rotor	5000	4x1x50ml	4420
5E-5	Swing Rotor	5000	4x1x100ml	4420
5E-6	Swing Rotor	4000	4x2x50ml	2820
5E-7	Swing Rotor	4000	4x4x15ml	2820
5E-8	Swing Rotor	4000	4x6x10/7ml vacuum tube	2580
5E-9	Swing Rotor	4000	4x6x5 ml vacuum tube adapter	2170